## Data Sheet



CellGenix<sup>®</sup> Recombinant Human FMS-like Tyrosine Kinase 3 Ligand (rh Flt-3L)

Preclinical Grade - Order No.: 1415-010 (10 μg), 1415-050 (50 μg)

### **Product Characteristics**

Source	E. coli
Description	Truncated human Flt-3L, accession # P49771, Thr27-Ala181 N-terminal Met and C-terminal 6xHis-tag Molecular mass 18.6 kDa
Formulation	Lyophilized from a 0.2 $\mu$ m-filtered solution containing 1.5 mM potassium phosphate, 8.1 mM sodium phosphate, 2.7 mM potassium chloride, and 137 mM sodium chloride, pH 7.5
Intended use	For preclinical ex vivo use. Not intended for therapeutic use.

### **Quality Parameters**

Activity	$\geq$ 0.5 x 10 <sup>6</sup> IU/mg calibrated against NIBSC #96/532 Measured in a cell proliferation assay using an Flt-3L-dependent cell line, OCI-AML5
Purity	$\geq$ 95 %, as determined by SDS-PAGE (under reducing and non-reducing conditions, visualized by Coomassie staining)
Endotoxin	< 1000 EU/mg, as determined by LAL gel clot test
Sterility	Sterility test of the vialed product (direct inoculation)
Mass per vial	1415-010: 10 μg, 1415-050: 50 μg
Animal-derived component-free	<b>ADCF Level 2:</b> The final product contains neither animal- nor human-derived materials. ADCF Level 2 cytokines are produced in our dedicated animal-free facility. No animal-derived components are used throughout the complete production process. All ADCF Level 2 cytokines are produced in <i>E. coli</i> .

# **S**CellGenix

### Shipment & Storage

ore in the
on. Store

Avoid repeated freeze/thaw cycles.

#### Handling Instructions

ReconstitutionRecommended in sterile water to a final concentration of 100 μg/ml (for 10 μg vials)<br/>or 250 μg/ml (for 50 μg vials).DilutionRecommended in CellGenix® serum-free media. For dilution with protein free<br/>medium, a carrier protein (0.1-1 % albumin or 1-10 % appropriate serum) has to be<br/>included. Failure to dilute product according to these instructions may result in loss<br/>of activity.

### **Quality Statement**

Final manufacturing steps and QC are performed in a GMP facility.

ZY-DS-1415d | 09/2017

HQ: CellGenix GmbH | Am Flughafen 16 | D - 79108 Freiburg | GERMANY | Phone +49 761 888 89-0 | Fax +49 761 888 89-830 US: CellGenix Inc. | One New Hampshire Avenue, Suite 125 | Portsmouth, NH 03801 | USA | Phone +1 603 373 0408 | Fax +1 603 373 8104 www.cellgenix.com | info@cellgenix.com